




中天医疗™  
Zenith Vascular

# Zenith Distal Access Catheter

NMPA No.20223031393

---

DEVELOPING INNOVATION  
TECHNOLOGIES TO OPTIMIZE  
VASCULAR TREATMENT





- **Navigation and Trackability**

Lubricious coating and smooth catheter transitions allow for easy access to MCA and beyond.

- Ultral-soft tip
- Large .070" Lumen

Serial Number	Product code	OD	ID	Maximum guidewire compatible	Working length
Distal Access Catheter 1.0 [NMPA NO.20203030818]	B1555115	5F	0.055"	≤0.035"	115cm
	B1555125	5F	0.055"	≤0.035"	125cm
	B1670105	6F	0.070"	≤0.035"	105cm
	B1670115	6F	0.070"	≤0.035"	115cm
	B1670125	6F	0.070"	≤0.035"	125cm
Distal Access Catheter 2.0 [NMPA NO.202223031393]	B3555115	5F	0.055"	≤0.035"	115cm
	B3555125	5F	0.055"	≤0.035"	125cm
	B3555133	5F	0.055"	≤0.035"	133cm
	B3670105	6F	0.070"	≤0.035"	105cm
	B3670115	6F	0.070"	≤0.035"	115cm
	B3670125	6F	0.070"	≤0.035"	125cm
	B3670133	6F	0.070"	≤0.035"	133cm



Zenith Vascular Scitech

TEL: 86-0512-6913873

E-mail: [Information@zenithmedtech.com](mailto:Information@zenithmedtech.com)

Website: [Zenithmedtech.com](http://Zenithmedtech.com)