





中天医疗™
Zenith Vascular

Zenith Microcatheter


NMPA No.20223030593

DEVELOPING INNOVATION
TECHNOLOGIES TO OPTIMIZE
VASCULAR TREATMENT





中天医疗™
Zenith Vascular



● 80cm hydrophilic coating

● Balanced shaft transitions and robust column strength minimize tip movement, providing reliable retrieval of intervention devices

● Flat braid design with lubricious PTFE liner maximizes lumen ID, providing smooth, low-friction delivery of intervention devices.

● Soft tip.



Peel-able sheath



Shape Needles

Product Code	ID [inch/mm]	Distal OD [mm]	Proximal OD [mm]	Working length [cm]
200017150	0.017" / 0.43	0.66	0.81	150
200017153	0.017" / 0.43	0.66	0.81	153
200017155	0.017" / 0.43	0.66	0.81	155
200021150	0.021" / 0.53	0.81	0.91	150
200021153	0.021" / 0.53	0.81	0.91	153
200021155	0.021" / 0.53	0.81	0.91	155
200027150	0.027" / 0.69	0.93	0.93	150
200027155	0.027" / 0.69	0.93	0.93	155



Zenith Vascular Scitech

TEL: 86-0512-6913873

E-mail: Information@zenithmedtech.com

Website: Zenithmedtech.com